















COLD & FLU GOT YOU DOWN?

The 2017-2018 flu season is reaching moderately severe levels in almost all 50 states. This season the flu vaccine is less effective against the predominate H3N2 strain. According to the CDC the most effective thing you can do to protect yourself, besides the flu shot, is wash your hands regularly with soap and water and use alcohol based hand sanitizer.

Protect your business from absenteeism and keep your employees healthy by ensuring they have the right products to prevent the spread of germs!



Cover your cough & sneeze!

Facial Tissue

ITEM	DESCRIPTION	BRAND	CASE PACK
GP47580	2PLY FACIAL TISSUE	ANGEL SOFT	96/54 SHT/CS
KC21340	2PLY FACIAL TISSUE	SURPASS	30/100 SHT/CS



Kills 99.9% of germs on hands!!

Hand Sanitizer

ITEM	DESCRIPTION	BRAND	CASE PACK
GOJ965212	TABLE TOP PUMP BOTTLE	PURELL	12/8 OZ/CS
GOJ192004	WHITE TOUCH-FREE DISPENSER	PURELL	1/EA
GOJ190502	FOAM HAND SANITIZER FOR DISPENSER	PURELL	2/1200 ML/CS



Break the chain of infection!!!

Surface Disinfectant

ITEM	DESCRIPTION	BRAND	CASE PACK
GOJ334206	PROFESSIONAL SURFACE SPRAY	PURELL	6/32 OZ/CS
CLO15948	LEMON FRESH DISINFECTING WIPES	CLOROX	6/75 CT/CS
CLO15949	FRESH SCENT DISINFECTING WIPES	CLOROX	6/75 CT/CS
CLO30830	HYDROGEN PEROXIDE DISINFECTANT WIPE	CLOROX	6/110 CT/CS



Clorox® Total 360™System

Superior Coverage with Trusted Clorox Solutions

Clorox® is proud to offer the first system that pairs an electrostatic sprayer with a portfolio of Clorox® disinfectants and sanitizers. Electrostatic technology enables superior coverage of Clorox® trusted solutions in hard to reach places — the side, underside and backside of surfaces. Using the Clorox® Total 360™ Electrostatic Sprayer to apply Clorox® Total 360™ Disinfectant Cleaner and Clorox® Anywhere® gives you a healthier facility.

CLX 360 device CLO60010
CLX 360 disinfectant 128oz CLO31650
CLX Anywhere sanitizer 128oz CLO31651





Clorox[®] **Total 360**[™] **System provides up to:**

TIME EFFICIENT

4X
FASTER*

SUPERIOR COVERAGE

18,000

65%

FT./HOUR LESS SOLUTION*